

Reusable EEG Cup Electrode Gold and Silver



Technomed understands the importance of a crystal clear signal, ease of use and electrode durability.

Technomed's Reusable EEG Cup Electrode Gold and Silver are designed to provide reliable signals and longtime ease of use.

These electrodes are designed for easy adjustment to patient contours to avoid electrode lift. The connector is optimized for easy headbox insertion and extraction.

User Benefits

- Clear and reliable signal
- Short cup-to-cable connection to reduce lift risk
- Optimal patient comfort with low profile disk, soft and flexible anti-tangle lead wires
- Optimized connector design for easy headbox insertion and extraction
- Long electrode life due to watertight design and Kevlar reinforced wires
- Available in two different cup sizes

Available variants

Materials	Cup Size	Cup Height	Wire Length	Packing	Part Number
Silver plated (Ag)	10 mm	3.1mm	1 m	Sets of 12, 2 x6 colors	TE/C12-434
Silver plated (Ag)	10 mm	3.1 mm	1.5 m	Sets of 12, 2 x6 colors	TE/C12-634
Silver plated (Ag)	10 mm	3.1 mm	2 m	Sets of 12, 2 x6 colors	TE/C12-834
Silver plated (Ag)	10 mm	3.1 mm	2.5 m	Sets of 12, 2 x6 colors	TE/C12-934
Silver plated (Ag)	6 mm	1.7 mm	1 m	Sets of 12, 2 x6 colors	TE/C14-434
Silver plated (Ag)	6 mm	1.7 mm	1.5 m	Sets of 12, 2 x6 colors	TE/C14-634
Silver plated (Ag)	6 mm	1.7 mm	2 m	Sets of 12, 2 x6 colors	TE/C14-834
Gold plated (Au)	10 mm	3.1 mm	1 m	Sets of 12, 2 x6 colors	TE/C32-434
Gold plated (Au)	10 mm	3.1 mm	1.5 m	Sets of 12, 2 x6 colors	TE/C32-634
Gold plated (Au)	10 mm	3.1 mm	2 m	Sets of 12, 2 x6 colors	TE/C32-834
Gold plated (Au)	10 mm	3.1 mm	2.5 m	Sets of 12, 2 x6 colors	TE/C32-934
Gold plated (Au)	6 mm	1.7 mm	1 m	Sets of 12, 2 x6 colors	TE/C34-434
Gold plated (Au)	6 mm	1.7 mm	1.5 m	Sets of 12, 2 x6 colors	TE/C34-634
Gold plated (Au)	6 mm	1.7 mm	2 m	Sets of 12, 2 x6 colors	TE/C34-834

2 x 6 Colors: red, green, white, yellow, blue, brown

Specifications

Cup	Ag or Au plating
Lead-free	√
Latex-free	√
BPA-free lead wire	√
DEHP-free lead wire	√
RoHS compliant	√
Lead wire design	Tinsel wire
Lead wire core	Kevlar and Silver Plated Copper
Lead wire insulation	TPU
Spring loaded connector	Gold plated Brass / PVC
Connector	DIN 42802

Cleaning and disinfection: see IFU for details



Taking you further. Step by step.



Your Technomed Business Partner:



Technomed

Amerikalaan 71
6199 AE | Maastricht-Airport | The Netherlands
T +31 43 408 68 68 | info@technomed.nl

Technomed USA

Suite 4990 Constellation Drive
White Bear Township | MN 55127 | United States
T +1 651 288 0890 | info@technomed-usa.com

www.technomed.nl